

MAP3K14 protein (STJ504444)
STJ504444**GENERAL INFORMATION**

Product Type	Proteins
Short Description	MAP3K14-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9020
Gene Symbol	MAP3K14
Uniprot ID	M3K14_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	