

Myoglobin protein (STJ504467)

STJ504467

GENERAL INFORMATION

Product Type	Proteins
Short Description	Myoglobin-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:1, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4151
Gene Symbol	MB
Uniprot ID	MYG_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	