

Nicotinic Acetylcholine Receptor alpha 7 protein (STJ504471)

STJ504471

GENERAL INFORMATION

Product Type	Proteins
Short Description	Nicotinic Acetylcholine Receptor alpha 7-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:250
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25302
Gene Symbol	Chrna7
Uniprot ID	ACHA7_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	