

nmt55/p54nrb protein (STJ504484)
STJ504484**GENERAL INFORMATION**

Product Type	Proteins
Short Description	nmt55/p54nrb-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4841
Gene Symbol	NONO
Uniprot ID	NONO_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	