

Niemann Pick C1 protein (STJ504487)

STJ504487

GENERAL INFORMATION

Product Type	Proteins
Short Description	Niemann Pick C1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4864
Gene Symbol	NPC1
Uniprot ID	NPC1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	