

OSR2 protein (STJ504499)
STJ504499**GENERAL INFORMATION**

Product Type	Proteins
Short Description	OSR2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:750
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	107587
Gene Symbol	Osr2
Uniprot ID	OSR2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	