

Nkx2.2 protein (STJ504501)

STJ504501

GENERAL INFORMATION

Product Type	Proteins
Short Description	Nkx2.2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4821
Gene Symbol	NKX2-2
Uniprot ID	NKX22_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 200-273 on human homeobox protein Nkx-2.2
Immunogen Region	200-273
Specificity	
Immunogen Sequence	