

CDKN2A/p14ARF protein (STJ504505)

STJ504505

GENERAL INFORMATION

Product Type	Proteins
Short Description	CDKN2A/p14ARF-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
Uniprot ID	ARF_HUMAN
Immunogen	Synthetic peptide corresponding to C-terminal amino acids 94-111 of Human p14ARF. Sequence: QRLGQQPLPRRPGHDDGQ
Immunogen Region	94-111
Specificity	
Immunogen	
Sequence	