

P2Y1 protein (STJ504507)

STJ504507

GENERAL INFORMATION

Product Type	Proteins
Short Description	P2Y1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25265
Gene Symbol	P2ry1
Uniprot ID	P2RY1_RAT
Immunogen	Synthetic peptide taken within amino acid region 360-410 on rat P2Y1.
Immunogen Region	360-410
Specificity	
Immunogen Sequence	