

Adenylate cyclase 1 protein (STJ504512)

STJ504512

GENERAL INFORMATION

Product Type	Proteins
Short Description	Adenylate cyclase 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	432530
Gene Symbol	Adcy1
Uniprot ID	ADCY1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1030-1060 on rat/mouse Adenylyl Cyclase 1 protein.
Immunogen Region	1030-1060
Specificity	
Immunogen Sequence	