

ADCY5 protein (STJ504516)

STJ504516

GENERAL INFORMATION

Product Type	Proteins
Short Description	ADCY5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64532
Gene Symbol	Adcy5
Uniprot ID	ADCY5_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on rat Particulate Adenylyl Cyclase 5 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	