

ADCY6 protein (STJ504517)

STJ504517

GENERAL INFORMATION

Product Type	Proteins
Short Description	ADCY6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25289
Gene Symbol	Adcy6
Uniprot ID	ADCY6_RAT
Immunogen	Synthetic peptide corresponding to N-terminal amino acids 13-27 of rat Particulate Adenylyl Cyclase 6 protein. Peptide was covalently modified and conjugated to KLH. Sequence: DERKTAWGERNGQKR
Immunogen Region	13-27
Specificity	
Immunogen Sequence	