

**ADCY8 protein (STJ504519)**  
STJ504519**GENERAL INFORMATION**

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Proteins   |
| <b>Short Description</b> | ADCY8-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition                 |
| <b>Host/Source</b>       |  |
| <b>Reactivity</b>        |  |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            |  |
| <b>Gene Symbol</b>        |  |
| <b>Uniprot ID</b>         |  |
| <b>Immunogen</b>          | Synthetic peptide corresponding to unique amino acid sequence on human Particulate Adenylyl Cyclase 8 protein. |
| <b>Immunogen Region</b>   |  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |