

alpha Actinin 4 protein (Phospho) (STJ504522)

STJ504522

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-alpha Actinin 4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81
Gene Symbol	ACTN4
Uniprot ID	ACTN4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-110 on human alpha-actinin-4 protien. Phosphorylated Serine-99 residue.
Immunogen Region	60-110
Specificity	
Immunogen	
Sequence	