

PALB2 protein (STJ504524)
STJ504524**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PALB2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79728
Gene Symbol	PALB2
Uniprot ID	PALB2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1, 130-1, 186 on human PALB2 protein.
Immunogen Region	1130-1186
Specificity	
Immunogen Sequence	