

Pannexin 1 protein (STJ504526)

STJ504526

GENERAL INFORMATION

Product Type	Proteins
Short Description	Pannexin 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24145
Gene Symbol	PANX1
Uniprot ID	PANX1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region from human Pannexin 1 protein
Immunogen Region	
Specificity	
Immunogen Sequence	