

Parkin protein (STJ504532)

STJ504532

GENERAL INFORMATION

Product Type	Proteins
Short Description	Parkin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5071
Gene Symbol	PRKN
Uniprot ID	PRKN_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 390-440 on human Parkin protein.
Immunogen Region	390-440
Specificity	
Immunogen Sequence	