

PAX1 protein (STJ504533)

STJ504533

GENERAL INFORMATION

Product Type	Proteins
Short Description	PAX1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5075
Gene Symbol	PAX1
Uniprot ID	PAX1_HUMAN
Immunogen	Synthetic peptide corresponding to unique sequence within amino acid region 210-260 on human PAX1 protein.
Immunogen Region	210-260
Specificity	
Immunogen Sequence	