

## Pax2 protein (STJ504534)

STJ504534

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Pax2-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5076</a>
<b>Gene Symbol</b>	<a href="#">PAX2</a>
<b>Uniprot ID</b>	<a href="#">PAX2_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 1-50, 255-305, and 310-360 on human PAX2 protein.
<b>Immunogen Region</b>	1-50 255-305 and 310-360
<b>Specificity</b>	
<b>Immunogen Sequence</b>	