

PAX4 protein (STJ504536)

STJ504536

GENERAL INFORMATION

Product Type	Proteins
Short Description	PAX4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	PAX4
Uniprot ID	PAX4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300, corresponding to unique amino acid sequence on PAX4 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	