

PAX5 protein (STJ504539)

STJ504539

GENERAL INFORMATION

Product Type	Proteins
Short Description	PAX5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5079
Gene Symbol	PAX5
Uniprot ID	PAX5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300 on human PAX5 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	