

PAX2 protein (STJ504574)
STJ504574**GENERAL INFORMATION**

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | PAX2-protein if for use in research applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: 1:500 |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 5076 |
| Gene Symbol | PAX2 |
| Uniprot ID | PAX2_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |