

Chimaerin beta 2 protein (Phospho) (STJ504583)

STJ504583

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Chimaerin beta 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1123
Gene Symbol	CHN1
Uniprot ID	CHIN_HUMAN
Immunogen	Phosphorylated synthetic peptide corresponding to unique amino acid sequence on beta 2 Chimaerin protein.
Immunogen Region	
Specificity	
Immunogen Sequence	