

NT5C2 protein (Phospho) (STJ504585)

STJ504585

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-NT5C2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22978
Gene Symbol	NT5C2
Uniprot ID	5NTC_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on cNII protein with phosphorylated Serine 7.
Immunogen Region	
Specificity	
Immunogen Sequence	