

PDE10A12 protein (STJ504589)

STJ504589

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE10A12-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10846
Gene Symbol	PDE10A
Uniprot ID	PDE10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-120 on PDE10A12 protein.
Immunogen Region	50-120
Specificity	
Immunogen Sequence	