

## PDE10A2 protein (STJ504591)

STJ504591

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Proteins   |
| <b>Short Description</b> | PDE10A2-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition                   |
| <b>Host/Source</b>       |  |
| <b>Reactivity</b>        |  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">10846</a>   |
| <b>Gene Symbol</b>        | <a href="#">PDE10A</a>  |
| <b>Uniprot ID</b>         | <a href="#">PDE10_HUMAN</a>   |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 1-50 on human PDE10A2 protein. |
| <b>Immunogen Region</b>   | 1-50  |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |