

PDE11A protein (STJ504595)

STJ504595

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE11A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50940
Gene Symbol	PDE11A
Uniprot ID	PDE11_HUMAN
Immunogen	Synthetic peptide conjugated to KLH, corresponding to C terminal amino acids 454-468 of human PDE11A.
Immunogen Region	454-468
Specificity	
Immunogen Sequence	