

PDE1A protein (STJ504602)

STJ504602

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE1A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to unique amino acid sequence within amino acid region 51-101 on PDE1A protein. Common to all PDE1A variants.
Immunogen Region	51-101
Specificity	
Immunogen Sequence	