

PDE1C protein (STJ504609)

STJ504609

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE1C-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5137
Gene Symbol	PDE1C
Uniprot ID	PDE1C_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on PDE1C gene.
Immunogen Region	
Specificity	
Immunogen Sequence	