

PDE4A1 protein (STJ504618)

STJ504618

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4A1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25638
Gene Symbol	Pde4a
Uniprot ID	PDE4A_RAT
Immunogen	Synthetic peptide specific for PDE4A1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	