

PDE4A5 protein (STJ504627)

STJ504627

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4A5-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:750-1:1, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5141
Gene Symbol	PDE4A
Uniprot ID	PDE4A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	