

PDE4A8 protein (STJ504628)

STJ504628

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4A8-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5141
Gene Symbol	PDE4A
Uniprot ID	PDE4A_HUMAN
Immunogen	Synthetic peptide specific only for PDE4A8 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	