

PDE4B3 protein (STJ504634)

STJ504634

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4B3-protein if for use in research applications.
Applications	ELISA, IF, IHC, IP, WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24626
Gene Symbol	Pde4b
Uniprot ID	PDE4B_RAT
Immunogen	Synthetic peptide selective only for PDE4B3 subtype.
Immunogen Region	
Specificity	
Immunogen Sequence	