

PDE4C protein (STJ504638)

STJ504638

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | PDE4C-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 5143 |
| Gene Symbol | PDE4C |
| Uniprot ID | PDE4C_HUMAN |
| Immunogen | Synthetic cyclic carboxy-terminal peptide common to all PDE4C variants. |
| Immunogen Region | C-Term |
| Specificity | |
| Immunogen Sequence | |