

PDE4C protein (STJ504639)
STJ504639**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PDE4C-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:2, 500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5143
Gene Symbol	PDE4C
Uniprot ID	PDE4C_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	