

PDE4D1 protein (STJ504640)

STJ504640

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4D1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide from unique region near N-terminus of PDE4D1 variant.
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	