

PDE4D2 protein (STJ504642)

STJ504642

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4D2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on PDE4D2 variant gene.
Immunogen Region	
Specificity	
Immunogen Sequence	