

PDE5A protein (STJ504659)

STJ504659

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE5A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8654
Gene Symbol	PDE5A
Uniprot ID	PDE5A_HUMAN
Immunogen	Recombinant protein common to all PDE5A variants corresponding to human PDE5A protein.
Immunogen Region	
Specificity	
Immunogen Sequence	