

PDE6 Beta protein (STJ504660)

STJ504660

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE6 Beta-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5158
Gene Symbol	PDE6B
Uniprot ID	PDE6B_HUMAN
Immunogen	Synthetic cyclic peptide common to all PDE6 Beta variants
Immunogen Region	
Specificity	
Immunogen Sequence	