

PDE7A protein (STJ504666)

STJ504666

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE7A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5150
Gene Symbol	PDE7A
Uniprot ID	PDE7A_HUMAN
Immunogen	Synthetic cyclic peptide common to PDE7A1 and PDE7A2 variants.
Immunogen Region	
Specificity	
Immunogen Sequence	