

PDE7B protein (STJ504671)

STJ504671

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE7B-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27115
Gene Symbol	PDE7B
Uniprot ID	PDE7B_HUMAN
Immunogen	Synthetic cyclic peptide common to all PDE7B variants.
Immunogen Region	
Specificity	
Immunogen Sequence	