

PDE9A protein (STJ504679)

STJ504679

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE9A-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: Apply 15-20 µl per lane
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5152
Gene Symbol	PDE9A
Uniprot ID	PDE9A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	