

Eph receptor A4 protein (Phospho) (STJ504687)

STJ504687

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Eph receptor A4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13838
Gene Symbol	Epha4
Uniprot ID	EPHA4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 570-630 on human EPHA R4 protein.
Immunogen Region	570-630
Specificity	
Immunogen Sequence	