

PERK protein (STJ504689)

STJ504689

GENERAL INFORMATION

Product Type	Proteins
Short Description	PERK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9451
Gene Symbol	EIF2AK3
Uniprot ID	E2AK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1000-1050 on human Eukaryotic translation initiation factor 2-alpha kinase 3 protein.
Immunogen Region	1000-1050
Specificity	
Immunogen Sequence	