

FGFR4 protein (Phospho) (STJ504691)

STJ504691

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-FGFR4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2264
Gene Symbol	FGFR4
Uniprot ID	FGFR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 530-580 on human Fibroblast growth factor receptor 4 variant.
Immunogen Region	530-580
Specificity	
Immunogen Sequence	