

PGAP3 protein (STJ504695)

STJ504695

GENERAL INFORMATION

Product Type	Proteins
Short Description	PGAP3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	93210
Gene Symbol	PGAP3
Uniprot ID	PGAP3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 10-60 on human PGAP3 protein.
Immunogen Region	10-60
Specificity	
Immunogen Sequence	