

## PGCA protein (STJ504698)

STJ504698

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	PGCA-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	4881
<b>Gene Symbol</b>	NPR1
<b>Uniprot ID</b>	ANPRA_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to C terminal amino acids 1042-1057 of Human PGCA. Sequence: VRTYWLLGERGCSTRG
<b>Immunogen Region</b>	1042-1057
<b>Specificity</b>	
<b>Immunogen Sequence</b>	