

PGCB protein (STJ504700)

STJ504700

GENERAL INFORMATION

Product Type	Proteins
Short Description	PGCB-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4882
Gene Symbol	NPR2
Uniprot ID	ANPRB_HUMAN
Immunogen	Synthetic C-terminal peptide corresponding to unique amino acid sequences on Particulate Guanylyl Cyclase B protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	