

Particulate Guanylate Cyclase C Peptide protein (STJ504702)

STJ504702

GENERAL INFORMATION

Product Type	Proteins
Short Description	Particulate Guanylate Cyclase C Peptide-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25711
Gene Symbol	Gucy2c
Uniprot ID	GUC2C_RAT
Immunogen	Synthetic C-terminal peptide corresponding to unique amino acid sequences on Particulate Guanylyl Cyclase C protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	