

PGCG protein (STJ504710)

STJ504710

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | PGCG-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | NA |
| Gene Symbol | Gucy2g |
| Uniprot ID | GUC2G_RAT |
| Immunogen | Synthetic C-terminal peptide corresponding to unique amino acid sequences on Particulate Guanylyl Cyclase G protein. |
| Immunogen Region | C-Term |
| Specificity | |
| Immunogen Sequence | |